


Product Datasheet

FBL Antibody (orb1248652)

Description	FBL Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: CRGKEDALVTKN	orb1248652 (0.1 ug/ml) staining of HEK29...
Target	FBL	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 35kDa band observed in lysates of cell line HEK293 (calculated MW of 33.8kDa according to NP_001427.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation 1 hour at room temperature. Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.	
Clonality	Polyclonal	
MW	Approx 34 kDa	
Uniprot ID	P22087	
NCBI	NP_001427.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 35kDa band observed in lysates of cell line HEK293 (calculated MW of 33.8kDa according to NP_001427.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation 1 hour at room temperature. Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.	
Expiration Date	12 months from date of receipt.	