

Product Datasheet

TAS2R4 Antibody (orb1248656)

Description

TAS2R4 Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

The immunogen for this antibody is: C-PELVTRNNTSFNIS

Target

TAS2R4

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Clonality

Polyclonal

Uniprot ID

[Q9NYW5](#)

NCBI

[NP_058640.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?