

Product Datasheet

USP6 Antibody (orb1248696)

Description

USP6 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications

Immunogen

The immunogen for this antibody is: C-HHIDL DVRTTLRNH

Target

USP6

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Preliminary experiments in Human Testis, Prostate and Placenta lysates gave no specific signal but low background at antibody concentration up to 1ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[P35125](#)

NCBI


[NP_004496.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Preliminary experiments in Human Testis, Prostate and Placenta lysates gave no specific signal but low background at antibody concentration up to 1ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

Expiration Date

12 months from date of receipt.



orb1248696
Immunofluorescence
analysis o...



orb1248696
Immunofluorescence
analysis o...



orb1248696 Flow
cytometric analysis
of p...