

Product Datasheet

cuff Antibody (orb1248711)

Description

cuff Antibody

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The immunogen for this antibody is: C-RRDENAHHVHNIEKT
Target	cuff
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Clonality	Polyclonal
Uniprot ID	Q9V629
NCBI	NP_610709.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Expiration Date	12 months from date of receipt.