







Product Datasheet

TNFRSF1A Antibody (orb1248726)

Description	TNFRSF1A Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1248726 (0.3 ug/ml) staining of Human...
Immunogen	The immunogen for this antibody is: C-SDHEIDRLELQNGR	
Target	TNFRSF1A	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248726 (0.5 ug/ml) staining of A549 ...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 52kDa band observed in Human Lung lysates and in lysates of cell lines A549, HepG2 and K562 (calculated MW of 50.5kDa according to NP_001056.1). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
MW	Approx 52 kDa	
Uniprot ID	P19438	
NCBI	NP_001056.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 52kDa band observed in Human Lung lysates and in lysates of cell lines A549, HepG2 and K562 (calculated MW of 50.5kDa according to NP_001056.1). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	