


Product Datasheet

MAPK9 Antibody (orb1248729)

Description	MAPK9 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1248729 (0.03 ug/ml) staining of HeLa...
Immunogen	The immunogen for this antibody is: ELVKGCVIFQGDH	
Target	MAPK9	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 48kDa band observed in lysates of cell line HeLa (calculated MW of 48.1kDa according to NP_002743.3). Recommended concentration: 0.03-0.1ug/ml.	
Clonality	Polyclonal	
MW	Approx 48 kDa	
Uniprot ID	P45984	
NCBI	NP_002743.3, NP_620707.1, NP_001128516.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 48kDa band observed in lysates of cell line HeLa (calculated MW of 48.1kDa according to NP_002743.3). Recommended concentration: 0.03-0.1ug/ml.	
Expiration Date	12 months from date of receipt.	