


Product Datasheet

GRN Antibody (orb1248740)

Description	GRN Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested	ELISA, WB	
Applications		orb1248740
Immunogen	The immunogen for this antibody is: PQDTVCDLIQSK	(1 ug/ml)
Target	GRN	staining of
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	A431 ly...
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 75kDa band observed in lysates of cell lines A431, A549 and HeLa (calculated MW of 63.5kDa according to NP_002078.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (He and Bateman, Cancer Res. 1999 Jul 1;59(13):3222-9. PMID: 10397269). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
MW	Approx 75 kDa	
Uniprot ID	P28799	
NCBI	NP_002078.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 75kDa band observed in lysates of cell lines A431, A549 and HeLa (calculated MW of 63.5kDa according to NP_002078.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (He and Bateman, Cancer Res. 1999 Jul 1;59(13):3222-9. PMID: 10397269). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	