

Product Datasheet

ABCC9 Antibody (orb1248747)

Description

ABCC9 Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

The immunogen for this antibody is: C-TSEYSINNTGKAD

Target

ABCC9

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave an approx 90kDa band in Human Testis lysates and 100kDa in Rat Testis lysates after 0.1ug/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 170kDa according to NP_005682.2. The 90kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Clonality

Polyclonal

MW

Approx 90 kDa

Uniprot ID

[O60706](#)

NCBI

[NP_064693.2](#), [NP_005682.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave an approx 90kDa band in Human Testis lysates and 100kDa in Rat Testis lysates after

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

peptide. We would appreciate any feedback