

Product Datasheet

APP Antibody (orb1248748)

Description	APP Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: DAEFRHDSGYEVHHQK-C	orb1248748 (1 ug/ml) staining of Human C...
Target	APP	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 75kDa band observed in Human Brain (Cerebral and Frontal cortex) lysates (calculated MW of 78.7kDa according to NP_958817.1 and 72.6 according to NP_001129601.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour.	
Clonality	Polyclonal	
MW	Approx 75 kDa	
Uniprot ID	P05067	
NCBI	NP_958816.1 , NP_001129601.1 , NP_000475.1 , NP_958817.1 , NP_001129602.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 75kDa band observed in Human Brain (Cerebral and Frontal cortex) lysates (calculated MW of 78.7kDa according to NP_958817.1 and 72.6 according to NP_001129601.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour.	
Expiration Date	12 months from date of receipt.	