

Product Datasheet

PCBP3 Antibody (orb1248774)

Description

PCBP3 Antibody

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The immunogen for this antibody is: PHSTLSTLSHHPQPQ-C
Target	PCBP3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Testis, Adipose, Lung and several cell line lysates gave no specific signal but low background (at antibody concentration up to 2ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Clonality	Polyclonal
Uniprot ID	P57721
NCBI	NP_065389.2 , NP_001123613.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Testis, Adipose, Lung and several cell line lysates gave no specific signal but low background (at antibody concentration up to 2ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Expiration Date	12 months from date of receipt.