



## Product Datasheet

### UBASH3A Antibody (orb1248862)

| Description         | UBASH3A Antibody                                                                                                                                                                                                                                                                                                                                                                       |  250kDa<br>150kDa<br>100kDa |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Species/Host        | Goat                                                                                                                                                                                                                                                                                                                                                                                   |  50kDa<br>37kDa             |
| Reactivity          | Human                                                                                                                                                                                                                                                                                                                                                                                  | 25kDa<br>20kDa                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                                                                                                                                                           | 15kDa                                                                                                          |
| Tested Applications | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                                                         | orb1248862<br>(1 ug/ml)<br>staining of<br>Human S...                                                           |
| Immunogen           | The immunogen for this antibody is: C-DSGIRISSVFASP                                                                                                                                                                                                                                                                                                                                    |                                                                                                                |
| Target              | UBASH3A                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                |
| Preservatives       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                            |                             |
| Form/Appearance     | Liquid                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                |
| Concentration       | 500 ug/mL                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                |
| Storage             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                                                                                                                             | orb1248862<br>(3.8 ug/ml)<br>staining of<br>paraf...                                                           |
| Note                | For research use only                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                |
| Application notes   | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 69.8kDa according to NP_001001895.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Breast shows perinuclear staining in the epithelial cells of the lobules. Recommended concentration, 3-5ug/ml. |                                                                                                                |
| Clonality           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                             |                                                                                                                |
| Uniprot ID          | <a href="#">P57075</a>                                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                |
| NCBI                | <a href="#">NP_001001895.1</a> , <a href="#">NP_061834.1</a>                                                                                                                                                                                                                                                                                                                           |                                                                                                                |
| Dilution Range      | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 69.8kDa according to NP_001001895.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Breast shows perinuclear staining in the epithelial cells of the lobules. Recommended concentration, 3-5ug/ml. |                                                                                                                |
| Expiration Date     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                |