


Product Datasheet

ARHGDIG Antibody (orb1248888)


Description

ARHGDIG Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: C-EAVPEYRAPGRK
Target	ARHGDIG
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 25kDa band observed in lysates of cell line HeLa (calculated MW of 25.1kDa according to NP_001167.2). Recommended concentration: 0.5-2ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Brain (Cortex). Recommended concentration: 3.75ug/ml.
Clonality	Polyclonal
Uniprot ID	Q99819
NCBI	NP_001167.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 25kDa band observed in lysates of cell line HeLa (calculated MW of 25.1kDa according to NP_001167.2). Recommended concentration: 0.5-2ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Brain (Cortex). Recommended concentration: 3.75ug/ml.
Expiration Date	12 months from date of receipt.



orb1248888
(0.5 ug/ml)
staining of
HeLa ...



orb1248888
(3.75 ug/ml)
staining of
para...