

Product Datasheet

FOXC1 Antibody (orb1248947)

Description

FOXC1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC

Applications

Immunogen

The immunogen for this antibody is: RTSGAFVYDCSKF

Target

FOXC1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the cytoplasm and nucleus of U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of HEK293 cells. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q12948](#)

NCBI



[NP_001444.2](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the cytoplasm and nucleus of U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of HEK293 cells. Recommended concentration: 10ug/ml.

Expiration Date


12 months from date of receipt.



orb1248947
Immunofluorescence
analysis o...



orb1248947 (3.75
ug/ml) staining of
para...



orb1248947 (3.75
ug/ml) staining of
para...