


Product Datasheet

AADAT Antibody (orb1248964)

Description	AADAT Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: C-KPEDAKNPQKNTPK
Target	AADAT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 45kDa band observed in Human Liver lysates (calculated MW of 47.4kDa according to NP_057312.1). Recommended concentration: 1-3ug/ml.
Clonality	Polyclonal
Uniprot ID	Q8N5Z0
NCBI	NP_057312.1, NP_872603.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 45kDa band observed in Human Liver lysates (calculated MW of 47.4kDa according to NP_057312.1). Recommended concentration: 1-3ug/ml.
Expiration Date	12 months from date of receipt.



orb1248964
(1 ug/ml)
staining of
Human L...