



Product Datasheet

NTN1 Antibody (orb1248973)

Description	NTN1 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	The immunogen for this antibody is: KFQQREKKGKCKKA
Target	NTN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000.Immunofluorescence: Strong expression of the protein seen in the membrane and cytoplasm of U2OS cells. Recommended concentration: 10ug/ml.
Clonality	Polyclonal
Uniprot ID	095631
NCBI	NP_004813.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000.Immunofluorescence: Strong expression of the protein seen in the membrane and cytoplasm of U2OS cells. Recommended concentration: 10ug/ml.
Expiration Date	12 months from date of receipt.



orb1248973 (3.75 ug/ml) staining of para...



orb1248973 Immunofluorescence analysis o...