

Product Datasheet

Dennd4c Antibody (orb1248985)

Description	Dennd4c Antibody	 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa
Reactivity	Mouse	 25kDa 20kDa
Conjugation	Unconjugated	 15kDa
Tested Applications	ELISA, WB	orb1248985 (0.5 ug/ml) staining of Mouse...
Immunogen	The immunogen for this antibody is: C-RFKPQSPYRAYKDRS	
Target	Dennd4c	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 210kDa band observed in Mouse Spleen lysates (calculated MW of 211kDa according to NP_908976.1). Recommended concentration: 0.5-2ug/ml.	
Clonality	Polyclonal	
Uniprot ID	A6H8H2	
NCBI	NP_001074483.1 , NP_908976.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 210kDa band observed in Mouse Spleen lysates (calculated MW of 211kDa according to NP_908976.1). Recommended concentration: 0.5-2ug/ml.	
Expiration Date	12 months from date of receipt.	