




Product Datasheet

ALDH18A1 Antibody (orb1249026)

Description	ALDH18A1 Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-SEHGSLKYLH	orb1249026 (1 ug/ml) staining of MCF7 ly...
Target	ALDH18A1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1249026 (0.3 ug/ml) staining of NIH3T...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 85kDa band observed in lysates of cell lines NIH3T3 and MCF7 (calculated MW of 87.3kDa according to NP_002851.2). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P54886	
NCBI	NP_002851.2 , NP_001017423.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 85kDa band observed in lysates of cell lines NIH3T3 and MCF7 (calculated MW of 87.3kDa according to NP_002851.2). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	