

Product Datasheet

41886 Antibody (orb1249028)

Description

41886 Antibody

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1249028
(0.5 ug/ml)
staining of
Human...

Immunogen

The immunogen for this antibody is: C-ELSKFVKDFSGN

Target

41886

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 52kDa band observed in Human Brain (Cerebellum) and Mouse Brain lysates (calculated MW of 53.0kDa according to NP_536341.1). Recommended concentration: 0.3-1ug/ml.

Clonality

Polyclonal

Uniprot ID
[O43236](#)
NCBI

[NP_536340.1](#), [NP_536341.1](#), [NP_004565.1](#),
[NP_001185642.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 52kDa band observed in Human Brain (Cerebellum) and Mouse Brain lysates (calculated MW of 53.0kDa according to NP_536341.1). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.

