


Product Datasheet

GRAMD3 Antibody (orb1249045)

Description			
Species/Host	Goat		
Reactivity	Human, Mouse		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB	orb1249045 (1 ug/ml) staining of Human T...	
Immunogen	The immunogen for this antibody is: C-QTPESNSRD		
Target	GRAMD3		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:4000.Western Blot:Approx 45kDa band observed in Human Brain (Temporal Cortex) and Mouse Brain lysates (calculated MW of 47.9kDa according to NP_076416.2). Recommended concentration: 1-3ug/ml.		
Clonality	Polyclonal		
Uniprot ID	Q96HH9		
NCBI	NP_001139793.1, NP_001139792.1, NP_001139794.1, NP_076416.2, NP_001139791.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000.Western Blot:Approx 45kDa band observed in Human Brain (Temporal Cortex) and Mouse Brain lysates (calculated MW of 47.9kDa according to NP_076416.2). Recommended concentration: 1-3ug/ml.		
Expiration Date	12 months from date of receipt.		