

Product Datasheet

RDH5 Antibody (orb1249095)

Description

RDH5 Antibody

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

orb1249095
(2 ug/ml)
staining of
Mouse E...

Immunogen

The immunogen for this antibody is: C-RDRQSLPASNAFV

Target

RDH5

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 35kDa band observed in Mouse Eye lysates (calculated MW of 34.8kDa according to Mouse NP_598767.1). Recommended concentration: 2-4ug/ml. Immunohistochemistry: In paraffin embedded Human Retina shows staining of the Ganglion Cell Layer. Recommended concentration, 5-10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q92781](#)

NCBI


[NP_002896.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 35kDa band observed in Mouse Eye lysates (calculated MW of 34.8kDa according to Mouse NP_598767.1). Recommended concentration: 2-4ug/ml. Immunohistochemistry: In paraffin embedded Human Retina shows staining of the Ganglion Cell Layer. Recommended concentration, 5-10ug/ml.

Expiration Date


12 months from date of receipt.



250kDa
150kDa
100kDa



50kDa
37kDa
25kDa
20kDa
15kDa



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

orb1249095
(5 ug/ml)
staining of
paraffi...