

Product Datasheet

VEGF antibody (PE) (orb124910)

Description

Rabbit polyclonal antibody against VEGF

Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human VEGF
Target	VEGFA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124910 is the PE conjugated form of orb11554
Isotype	IgG
Clonality	Polyclonal
MW	23kDa
Uniprot ID	P15692
Dilution Range	Flow-Cyt=,IF=1:100-500
Expiration Date	12 months from date of receipt.