

## Product Datasheet

### HOXA9 Antibody (orb1249107)

**Description**

HOXA9 Antibody

**Species/Host**

Goat

**Reactivity**

Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

orb1249107  
(1 ug/ml)  
staining of  
Mouse a...

**Immunogen**

The immunogen for this antibody is: C-QPPRQAATLAHPD

**Target**

HOXA9

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Form/Appearance**

Liquid

**Concentration**

500 ug/mL

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 30 to 33kDa bands observed in Mouse and Rat Spinal Cord lysates (calculated MW of 29.9kDa according to Mouse NP\_034586.1 and 25.8kDa according to Rat NP\_001102703.2). Recommended concentration: 1-3ug/ml.

**Clonality**

Polyclonal

**Uniprot ID**

[P31269](#)

**NCBI**

[NP\\_689952.1](#)

**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 30 to 33kDa bands observed in Mouse and Rat Spinal Cord lysates (calculated MW of 29.9kDa according to Mouse NP\_034586.1 and 25.8kDa according to Rat NP\_001102703.2). Recommended concentration: 1-3ug/ml.

**Expiration Date**

12 months from date of receipt.