

Product Datasheet

POU3F3 Antibody (orb1249115)

Description

POU3F3 Antibody

250kDa
150kDa
100kDa

Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1249115
(1 ug/ml)
staining of
Mouse S...

Immunogen

The immunogen for this antibody is: C-HMLSHAHQWVTAL

Target

POU3F3

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 48kDa band observed in Mouse Spinal Cord lysates (calculated MW of 50.3kDa according to Human NP_006227.1 and 50.2kDa according to Mouse NP_032926.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour.

Clonality

Polyclonal

Uniprot ID

[P20264](#)

NCBI

[NP_006227.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 48kDa band observed in Mouse Spinal Cord lysates (calculated MW of 50.3kDa according to Human NP_006227.1 and 50.2kDa according to Mouse NP_032926.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour.

Expiration Date

12 months from date of receipt.