

Product Datasheet

ARMS2 Antibody (orb1249147)

Description	ARMS2 Antibody	 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa
Reactivity	Human	 25kDa 20kDa
Conjugation	Unconjugated	 15kDa
Tested Applications	ELISA, IHC	orb1249147 (4 ug/ml) staining of paraffi...
Immunogen	The immunogen for this antibody is: C-EGASDKQRSK	
Target	ARMS2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	 250kDa 150kDa 100kDa 75kDa 50kDa 37kDa
Form/Appearance	Liquid	25kDa 20kDa
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1249147 Negative Control showing stai...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 11.4kDa according to NP_001093137.1. Therefore we cannot recommend an optimal concentration. Immunohistochemistry: In paraffin embedded Human Retina shows cytoplasm staining throughout. Recommended concentration: 2-4ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P0C7Q2	
NCBI	NP_001093137.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 11.4kDa according to NP_001093137.1. Therefore we cannot recommend an optimal concentration. Immunohistochemistry: In paraffin embedded Human Retina shows cytoplasm staining throughout. Recommended concentration: 2-4ug/ml.	
Expiration Date	12 months from date of receipt.	