


Product Datasheet

RSL1D1 Antibody (orb1249183)

Description	RSL1D1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-PTARKQLDKEQVRK	orb1249183 (1 ug/ml) staining of K562 ly...
Target	RSL1D1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 75kDa band observed in lysates of cell line K562 (calculated MW of 55.0kDa according to NP_056474.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O76021	
NCBI	NP_056474.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 75kDa band observed in lysates of cell line K562 (calculated MW of 55.0kDa according to NP_056474.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	