

Product Datasheet

ACPP Antibody (orb1249203)

Description

ACPP Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

The immunogen for this antibody is: C-YGIHKQKEKSRLQ

Target

ACPP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 44.6kDa according to NP_001090.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Patel et al, Clin Chem. 1986 Oct;32(10):1832-5. PMID: 3530530). Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows cytoplasm staining in some secretory cells of the gland. Recommended concentration, 3-5ug/ml.

Clonality

Polyclonal

Uniprot ID

[P15309](#)

NCBI



[NP_001090.2](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 44.6kDa according to NP_001090.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Patel et al, Clin Chem. 1986 Oct;32(10):1832-5. PMID: 3530530). Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows cytoplasm staining in some secretory cells of the gland. Recommended concentration, 3-5ug/ml.

Expiration Date

12 months from date of receipt.



orb1249203
(0.3 ug/ml)
staining of
Human...



orb1249203
(3.8 ug/ml)
staining of
paraf...