


Product Datasheet

ACPP Antibody (orb1249269)

Description	ACPP Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-NHMKRATQIPSYKK	
Target	ACPP	orb1249269 (0.01 ug/ml) staining of Huma...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 48.3kDa according to NP_001127666.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows reticulate cytoplasm staining, more intens towards the lumen of the gland. Recommended concentration, 3-5ug/ml.	orb1249269 (3.8 ug/ml) staining of paraf...
Clonality	Polyclonal	
Uniprot ID	P15309	
NCBI	NP_001090.2, NP_001127666.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 48.3kDa according to NP_001127666.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows reticulate cytoplasm staining, more intens towards the lumen of the gland. Recommended concentration, 3-5ug/ml.	
Expiration Date	12 months from date of receipt.	