

Product Datasheet

C22orf28 Antibody (orb1249305)

Description

C22orf28 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

The immunogen for this antibody is: C-QADPNKVSARAKKR

Target

C22orf28

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 55kDa band observed in lysates of cell lines HeLa, Jurkat and K562 (calculated MW of 55.2kDa according to NP_055121.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis shows nuclear and vicinity staining in spermatogonia. Recommended concentration, 5-10ug/ml. Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.

Clonality



Polyclonal

Uniprot ID
[Q9Y3I0](#)
NCBI
[NP_055121.1](#)
Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 55kDa band observed in lysates of cell lines HeLa, Jurkat and K562 (calculated MW of 55.2kDa according to NP_055121.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis shows nuclear and vicinity staining in spermatogonia. Recommended concentration, 5-10ug/ml. Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.

Expiration Date

12 months from date of receipt.



orb1249305
(0.3 ug/ml)
staining of
HeLa ...



orb1249305
(5 ug/ml)
staining of
paraffi...