

Product Datasheet

TTLL4 Antibody (orb1249355)

Description

TTLL4 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

orb1249355
(5 ug/ml)
staining of
paraffi...

Immunogen

The immunogen for this antibody is: C-DSQDPTKEIRFTE

Target

TTLL4

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Immunohistochemistry: In paraffin embedded Human Lung shows cytoplasm staining in the respiratory epithelium of the bronchus. Recommended concentration, 5-10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q14679](#)

NCBI

[NP_055455.3](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Immunohistochemistry: In paraffin embedded Human Lung shows cytoplasm staining in the respiratory epithelium of the bronchus. Recommended concentration, 5-10ug/ml.

Expiration Date

12 months from date of receipt.

