

Product Datasheet

GHITM Antibody (orb1249360)

Description

GHITM Antibody



Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1249360
(0.01 ug/ml)
staining of
Huma...

Immunogen

The immunogen for this antibody is: C-PSREYATKTRIGIR

Target

GHITM

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 35kDa band observed in lysates of cell line HeLa and of Human Brain (Cerebellum), Duodenum and Heart (calculated MW of 37.2kDa according to Human NP_055209.2). Recommended concentration: 0.01-0.1ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q9H3K2](#)

NCBI

[NP_055209.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 35kDa band observed in lysates of cell line HeLa and of Human Brain (Cerebellum), Duodenum and Heart (calculated MW of 37.2kDa according to Human NP_055209.2). Recommended concentration: 0.01-0.1ug/ml.

Expiration Date

12 months from date of receipt.