


Product Datasheet

SPG20 Antibody (orb1249374)


Description

SPG20 Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The immunogen for this antibody is: C-EASGTDVKQLDQGNK
Target	SPG20
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 72.8kDa according to NP_055902.1). Recommended concentration: 1-3ug/ml. Immunofluorescence: HeLa show diminished signals after knock down of expression by introduction of siRNA.
Clonality	Polyclonal
Uniprot ID	Q8N0X7
NCBI	NP_055902.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 72.8kDa according to NP_055902.1). Recommended concentration: 1-3ug/ml. Immunofluorescence: HeLa show diminished signals after knock down of expression by introduction of siRNA.
Expiration Date	12 months from date of receipt.



orb1249374
(2 ug/ml)
staining of
HeLa ly...



orb1249374
(2 ug/ml)
staining of
methano...