

Product Datasheet

FBN1 Antibody (orb1249415)

Description
FBN1 Antibody
Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications
Immunogen

The immunogen for this antibody is: C-DASNIEDQSETEAN

Target

FBN1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes


Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Preliminary experiments in Human Brain (Cerebellum), Heart, Kidney, Skeletal Muscle and Skin lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunofluorescence: Strong expression of the protein seen in the cytoplasm of U251 and Jurkat cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of Jurkat cells. Recommended concentration: 10ug/ml.

Clonality


Polyclonal

Uniprot ID
P35555
NCBI
NP_000129.3
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Preliminary experiments in Human Brain (Cerebellum), Heart, Kidney, Skeletal Muscle and Skin lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunofluorescence: Strong expression of the protein seen in the cytoplasm of U251 and Jurkat cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of Jurkat cells. Recommended concentration: 10ug/ml.



orb1249415
Immunofluorescence
analysis o...



orb1249415
Immunofluorescence
analysis o...



orb1249415 Flow
cytometric analysis
of p...