


Product Datasheet

SUMF1 Antibody (orb1249460)

Description		SUMF1 Antibody			
Species/Host	Goat				
Reactivity	Human, Mouse				
Conjugation	Unconjugated				
Tested Applications	ELISA, WB				
Immunogen	The immunogen for this antibody is: C-ETLNPKGPPSGKDR	orb1249460 (1 ug/ml) staining of Mouse P...			
Target	SUMF1				
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.				
Form/Appearance	Liquid				
Concentration	500 ug/mL				
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.				
Note	For research use only				
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 37-40kDa band observed in Mouse Pancreas, Mouse Eye and Human Kidney lysates (calculated MW of 40.6kDa according to Human NP_877437.2 and 40.6kDa according to Mouse NP_666049.2). Recommended concentration: 0.5-2ug/ml.				
Clonality	Polyclonal				
Uniprot ID	Q8NBK3				
NCBI	NP_877437.2				
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 37-40kDa band observed in Mouse Pancreas, Mouse Eye and Human Kidney lysates (calculated MW of 40.6kDa according to Human NP_877437.2 and 40.6kDa according to Mouse NP_666049.2). Recommended concentration: 0.5-2ug/ml.				
Expiration Date	12 months from date of receipt.				