


Product Datasheet

UBE2F Antibody (orb1249493)

Description	UBE2F Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1249493 (0.05 ug/ml) staining of Mous...
Immunogen	The immunogen for this antibody is: C-HLRDKEDFRNKVDD	
Target	UBE2F	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	In paraffin embedded Human Placenta show...
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 20kDa band observed in Mouse Liver lysates (calculated MW of 21.1kDa according to Human NP_542409.1 and to Mouse NP_080730.1). Recommended concentration: 0.05-0.2ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows membrane staining in trophoblasts-. Recommended concentration, 2.5-3.8ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q969M7	
NCBI	NP_542409.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 20kDa band observed in Mouse Liver lysates (calculated MW of 21.1kDa according to Human NP_542409.1 and to Mouse NP_080730.1). Recommended concentration: 0.05-0.2ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows membrane staining in trophoblasts-. Recommended concentration, 2.5-3.8ug/ml.	
Expiration Date	12 months from date of receipt.	