



Product Datasheet

Caspase 9 antibody (orb12495)


Description	Rabbit polyclonal antibody to Caspase 9
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Caspase-9(NP_001220.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47KDa/35KDa/37KDa
Uniprot ID	P55211
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of human colon tissue.