

## Product Datasheet

### JUNB Antibody (orb1249619)

## Description

JUNB Antibody

### Species/Host

Goat

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

The immunogen for this antibody is: C-TLKAENAGLSSTAG

### Target

JUNB

### Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

### Form/Appearance

Liquid

### Concentration

500 ug/mL

### Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

### Note

For research use only

### Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 35kDa band observed in Human Breast Cancer lysates (calculated MW of 35.9kDa according to NP\_002220.1). Recommended concentration: 0.3-1ug/ml.

### Clonality

Polyclonal

### Uniprot ID

[P17275](#)

### NCBI


[NP\\_002220.1](#)

### Dilution Range

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 35kDa band observed in Human Breast Cancer lysates (calculated MW of 35.9kDa according to NP\_002220.1). Recommended concentration: 0.3-1ug/ml.

### Expiration Date

12 months from date of receipt.



orb1249619  
(0.3 ug/ml)  
staining of  
Human...