


Product Datasheet

STK17A Antibody (orb1249714)

Description	STK17A Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, RNA-I, WB
Immunogen	STK17A (AAH47696, 301 a.a. ~ 414 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	STK17A
Preservatives	PBS, pH 7.2
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Store DRAK1 antibody at -20 °C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Application notes	DRAK1 antibody can be used in Western Blot, and immunohistochemistry starting at 5 µg/mL.
Isotype	IgG2a,k
Clonality	Monoclonal
Uniprot ID	Q9UEE5
NCBI	Q9UEE5
Dilution Range	DRAK1 antibody can be used in Western Blot, and immunohistochemistry starting at 5 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry of human placenta t...