


Product Datasheet

GSK3B Antibody (orb1249759)

Description	GSK3B Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	GSK3B antibody was raised against phospho amino acids 4 - 12 of GSK3B (Human).
Target	GSK3B
Preservatives	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Short term: 4°C. Long term: Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	GSK3B antibody can be used in ELISA starting at 1:15000 - 1:60000, Western Blot starting at 1:500 - 1:2000, and immunohistochemistry starting at 5 µg/mL.
Clonality	Polyclonal
Uniprot ID	P49841
NCBI	P49841
Dilution Range	GSK3B antibody can be used in ELISA starting at 1:15000 - 1:60000, Western Blot starting at 1:500 - 1:2000, and immunohistochemistry starting at 5 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of GSK3B i...