

Product Datasheet

LIN28B Antibody (orb1249839)

Description

LIN28B Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	KLH-conjugated synthetic peptide between amino acids 1-30 from N-terminal region of human LIN28B.
Target	LIN28B
Preservatives	PBS, 0.09% sodium azide
Concentration	batch dependent
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	LIN28B Antibody can be used for ELISA starting at 1:1000, Flow Cytometry starting at 1:10, Immunohistochemistry in paraffin embedded tissues starting at 10µg/mL, and Western Blotting starting at 1:50 - 1:100.
Clonality	Polyclonal
Uniprot ID	Q6ZN17
NCBI	Q6ZN17
Dilution Range	LIN28B Antibody can be used for ELISA starting at 1:1000, Flow Cytometry starting at 1:10, Immunohistochemistry in paraffin embedded tissues starting at 10µg/mL, and Western Blotting starting at 1:50 - 1:100.
Expiration Date	12 months from date of receipt.