

Product Datasheet

CD19 Antibody (orb1249843)

Description	CD19 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	KLH-conjugated synthetic peptide between amino acids 393-421 from C-terminal region of human CD19.
Target	CD19
Preservatives	PBS, 0.09% sodium azide
Concentration	batch dependent
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	CEP57 Antibody can be used for ELISA starting at 1:40000, and Immunohistochemistry in paraffin embeded tissues starting at 5 µg/mL,
Clonality	Polyclonal
Uniprot ID	P15391
NCBI	P15391
Dilution Range	CEP57 Antibody can be used for ELISA starting at 1:40000, and Immunohistochemistry in paraffin embeded tissues starting at 5 µg/mL,
Expiration Date	12 months from date of receipt.