

Product Datasheet

ACE2 Antibody (orb1249850)

Description	ACE2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	KLH-conjugated synthetic peptide between amino acids 140-172 from internal region of human ACE2 (SARS Receptor).
Target	ACE2
Preservatives	PBS, 0.09% sodium azide
Concentration	batch dependent
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	Flo (1:10 - 1:50), IHC-P (5 ug/ml), WB (1:1000)
Clonality	Polyclonal
Uniprot ID	Q9BYF1
NCBI	Q9BYF1
Dilution Range	Flo (1:10 - 1:50), IHC-P (5 ug/ml), WB (1:1000)
Expiration Date	12 months from date of receipt.