

Product Datasheet

MCL1 Antibody (orb1249876)

Description	MCL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide corresponding to aa 121-139 of human MCL-1.
Target	MCL1
Preservatives	10 mM PBS, pH 7.4, rAlbumin, sodium azide.
Concentration	batch dependent
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	MCL1 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:25 - 1:100.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07820
NCBI	Q07820
Dilution Range	MCL1 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:25 - 1:100.
Expiration Date	12 months from date of receipt.