

Product Datasheet

F9 Antibody (orb1249889)

Description	F9 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Full length native protein (purified) (Human).
Target	F9
Preservatives	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	F9 Monoclonal antibody may be stored at Store at -70°C
Note	For research use only
Application notes	F9 Monoclonal Antibody can be used for immunohistochemistry (tested on paraffin embedded tissues) starting at 10 µg/mL.
Isotype	IgG1,k
Clonality	Monoclonal
Uniprot ID	P00740
NCBI	P00740
Dilution Range	F9 Monoclonal Antibody can be used for immunohistochemistry (tested on paraffin embedded tissues) starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.