

Product Datasheet

NAPEPLD Antibody (orb1249898)

Description	NAPEPLD Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human NAPE-PLD amino acids 6-20
Target	NAPEPLD
Preservatives	TBS pH 7.4, containing 50% glycerol, 0.1% rAlbumin, and 0.02% sodium azide
Concentration	batch dependent
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	NAPEPLD / NAPE-PLD Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 5 µg/mL, and Western Blotting starting at 1:200.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IQ20
NCBI	Q6IQ20
Dilution Range	NAPEPLD / NAPE-PLD Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 5 µg/mL, and Western Blotting starting at 1:200.
Expiration Date	12 months from date of receipt.