

## Product Datasheet

### COL5A1 Antibody [Out of stock] (orb1249901)

**Description**

COL5A1 Antibody [Out of stock]

**Species/Host**

Rabbit

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IHC, IP, WB

**Target**

COL5A1

**Preservatives**

0.125 M sodium borate, 0.075 M sodium chloride, 0.005 M EDTA;, pH 8.0, 0.01% sodium azide.

**Concentration**

1 mg/ml

**Storage**

Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.

**Note**

For research use only

**Application notes**

Collagen V Antibody can be used for ELISA starting at 1:20000 - 1:100000, Immunohistochemistry in paraffin embeded tissues starting at 5 µg/mL, Immunoprecipitation and Western Blotting starting at 1:2000 - 1:10000.

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID**
**P20908**
**NCBI**
**P20908**
**Dilution Range**

Collagen V Antibody can be used for ELISA starting at 1:20000 - 1:100000, Immunohistochemistry in paraffin embeded tissues starting at 5 µg/mL, Immunoprecipitation and Western Blotting starting at 1:2000 - 1:10000.

**Expiration Date**

12 months from date of receipt.