


Product Datasheet

IL10 Antibody (orb1249903)

Description	IL10 Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant mouse IL-10
Target	IL10
Preservatives	PBS, 0.09% sodium azide. Sourced from tissus culture supernatant
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store IL-10 antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	ELISAFlo (1:10 - 1:100)IHC-FrIHC-P (10 ug/ml)WB (1:1000)
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P22301
NCBI	P22301
Dilution Range	ELISAFlo (1:10 - 1:100)IHC-FrIHC-P (10 ug/ml)WB (1:1000)
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of mouse small inte...