

Product Datasheet

CTNND1 Antibody (orb1249932)

Description	CTNND1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide corresponding to C-Terminal of human p120.
Target	CTNND1
Preservatives	10 mM PBS, pH 7.4, rAlbumin, sodium azide.
Concentration	batch dependent
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	CTNND1 / p120CTN Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:15.
Clonality	Polyclonal
Uniprot ID	O60716
NCBI	O60716
Dilution Range	CTNND1 / p120CTN Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:15.
Expiration Date	12 months from date of receipt.