

Product Datasheet

NRG1 Antibody (orb1249962)

Description	NRG1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide from human NRG1 isoform-10.
Target	NRG1
Preservatives	PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol
Concentration	1 mg/ml
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	ELISA (1:10000)IHC-P (1:200)WB (1:500 - 1:1000)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q02297
NCBI	Q02297
Dilution Range	ELISA (1:10000)IHC-P (1:200)WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.