

Product Datasheet

UBE2R2 antibody (orb125001)

Description

Goat polyclonal antibody to UBE2R2

Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IEM, IHC

Applications
Target

UBC3B

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: No signal obtained yet but low background observed in Human Testis and U937 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?. IHC: In paraffin embedded Human Testis shows textured nuclear staining of the primary spermatocytes . Recommended concentration, 2-4µg/ml.. Notes: Immunoelectron microscopy: This product stained chromatoid bodies in Rat spermatids.

Protein Sequence

AQQQMTSSQKALM

Clonality

Polyclonal

MW

27

Entrez


54926

Dilution Range

ELISA: 1:4000, WB: 1 µg/ml, IHC-P: 2-4µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Testis