


Product Datasheet

CD69 Antibody (orb1250063)

Description	CD69 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Recombinant protein
Target	CD69
Preservatives	0.03% sodium azide
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Short term: 4°C. Long term: Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	ELISAFLOW (1:200 - 1:400)IHC-P (1:200)WB (1:500 - 1:2000)
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q07108
NCBI	Q07108
Dilution Range	ELISAFLOW (1:200 - 1:400)IHC-P (1:200)WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.



Immunohistochemistry of human small intestine tissue.